

Smartez™ Infusion Pumps - Soft Shell

Short Infusion Duration

Model Number (REF)	Nominal Volume	Nominal Flow Rate	Nominal Time
481010	100 ml	200 ml/h	30 min
481020	250 ml	500 ml/h	30 min
481030	50 ml	50 ml/h	60 min
481040	100 ml	100 ml/h	60 min
481050	250 ml	250 ml/h	60 min
481060	250 ml	175 ml/h	90 min
481070	400 ml	200 ml/h	120 min
481080	500 ml	250 ml/h	120 min
481090	100 ml	50 ml/h	120 min
481100	400 ml	100 ml/h	240 min
481110	200 ml	100 ml/h	120 min
481120	250 ml	100 ml/h	150 min
481130	250 ml	125 ml/h	120 min
481140	200 ml	200 ml/h	60 min
481150	250 ml	50 ml/h	300 min
481160	120 ml	40 ml/h	180 min
481190	50 ml	100 ml/h	30 min
481200	250 ml	20 ml/h	12.5 hrs
481210	250 ml	21 ml/h	12 hrs
481220	120 ml	30 ml/h	4 hrs
481230	270 ml	175 ml/h	90 min
481250	200 ml	25 ml/h	8 hrs

Smartez™ Infusion Pumps - Soft Shell

Long Infusion Duration

Model Number (REF)	Nominal Volume	Nominal Flow Rate	Nominal Time
480010	60 ml	5 ml/h	12 hrs
480020	80 ml	5 ml/h	16 hrs
480030	125 ml	5 ml/h	25 hrs
480040	270 ml	10 ml/h	27 hrs
480050	60 ml	2 ml/h	30 hrs
480060	120 ml	4 ml/h	30 hrs
480070	400 ml	10 ml/h	40 hrs
480080	100 ml	2 ml/h	50 hrs
480090	270 ml	5 ml/h	54 hrs
480100	120 ml	2 ml/h	60 hrs
480110	400 ml	5 ml/h	80 hrs
480120	100 ml	1.5 ml/h	67 hrs
480130	270 ml	4 ml/h	68 hrs
480140	400 ml	4 ml/h	100 hrs
480150	65 ml	0.5 ml/h	130 hrs
480160	270 ml	2 ml/h	135 hrs
480170	300 ml	2 ml/h	150 hrs
480180	100 ml	0.5 ml/h	200 hrs
480190	270 ml	1 ml/h	270 hrs
480200	240 ml	2 ml/h	120 hrs
480210	100 ml	8 ml/h	12.5 hrs
480220	100 ml	4 ml/h	25 hrs
480230	270 ml	1.5 ml/h	180 hrs

Smartez™ Infusion Pumps - Soft Shell

Long Infusion Duration

Model Number (REF)	Nominal Volume	Nominal Flow Rate	Nominal Time
480250	120 ml	10 ml/h	12 hrs
480260	240 ml	10 ml/h	24 hrs
480270	250 ml	5 ml/h	50 hrs
480280	60 ml	2.5 ml/h	24 hrs
480290	250 ml	1.5 ml/h	167 hrs
480300	85 ml	0.5 ml/h	170 hrs
480310	100 ml	1 ml/h	100 hrs
480320	270 ml	3.75 ml/h	72 hrs
480330	270 ml	2.25 ml/h	120 hrs
480340	60 ml	1.5 ml/h	40 hrs
480350	300 ml	5 ml/h	60 hrs
480360	300 ml	10 ml/h	30 hrs
480370	250 ml	8 ml/h	31 hrs
480380	120 ml	1 ml/h	120 hrs
480390	150 ml	3 ml/h	50 hrs
480410	150 ml	6 ml/h	25 hrs
480420	275 ml	2.5 ml/h	110 hrs
480430	300 ml	7 ml/h	43 hrs
480440	300 ml	4 ml/h	75 hrs

Smartez™ Infusion Pumps - Soft Shell

Long Infusion Duration

Model Number (REF)	Nominal Volume	Nominal Flow Rate	Nominal Time
480450	300 ml	6 ml/h	50 hrs
480460	60 ml	0.5 ml/h	120 hrs
480470	60 ml	1 ml/h	60 hrs
480480	275 ml	0.5 ml/h	550 hrs
480490	275 ml	7 ml/h	43 hrs
480500	275 ml	1 ml/h	275 hrs
480510	275 ml	1.5 ml/h	183 hrs
480520	275 ml	2 ml/h	138 hrs
480530	275 ml	6.5 ml/h	42 hrs
480540	275 ml	3 ml/h	92 hrs
480550	275 ml	3.5 ml/h	79 hrs
480560	275 ml	4 ml/h	69 hrs
480570	275 ml	4.5 ml/h	61 hrs
480580	275 ml	5 ml/h	55 hrs
480590	275 ml	5.5 ml/h	50 hrs
480600	275 ml	6 ml/h	46 hrs
480610	100 ml	5 ml/h	20 hrs
480620	240 ml	5 ml/h	48 hrs
480630	84 ml	0.5 ml/h	168 hrs
480640	96 ml	2 ml/h	48 hrs

Smartez™ Infusion Pumps - Soft Shell

Long Infusion Duration

Model Number (REF)	Nominal Volume	Nominal Flow Rate	Nominal Time
480650	252 ml	1.5 ml/h	168 hrs
480660	48 ml	2 ml/h	24 hrs
480670	200 ml	2 ml/h	100 hrs
481170	500 ml	20 ml/h	25 hrs
481180	500 ml	40 ml/h	12.5 hrs